

AMA Adaptive Devices

Ensuring accessibility



We understand that some members have limitations, but we always find ways to help everyone. Using adaptive technology devices is important for making sure everyone can access and use our services. There are different ways to activate devices for calling for help, like pillow push pads or sip-and-puff devices. These can be customized for your clients. Our Adaptive Device solutions are only limited by what we can imagine, so we need your input to come up with the best solutions for everyone.

HERE ARE SOME FAVORITES:



Pillow Switches

The Pillow Switch is a wired assistive technology device that can be conveniently activated by a client's head or cheek, and easily attached to a pillow or wheelchair cushion with Velcro and a safety pin. It provides visual and auditory feedback when activated, making it ideal for users with limited physical movement or dexterity.



What is Sip and Puff?

The Sip and Puff (SNP) device allows users to control their AMA Personal Response System by sipping or puffing through a tube, straw, or wand using air pressure signals. The 36-inch Easy Flex Sip and Puff Switch is designed for individuals with limited mobility or physical challenges for easy mouth operation.



Bright Switch for the Visually Impaired

A great choice for visually impaired users, this switch's bright blinking light makes it easy to locate. When the switch is activated, its textured yellow light shines, music plays and the switch vibrates.

Is this an extra cost?

NO! AMA Adaptive Devices are complimentary with your client's PERS authorization.

